

Guardian[®] Class I Biological Safety Cabinets



MSC 1200 C1

**BRITISH
MANUFACTURED**



**FIVE YEAR
WARRANTY**



**MAXIMUM
OPERATOR
SAFETY**



**HEALTH
PROTECTION
AGENCY TESTED**



The Guardian[®] Class I Biological Safety Cabinet sets new standards for operator and environmental protection. Designed to incorporate the latest features to BSEN 12469:2000 standards.

Standard Features include:

- Low Energy <120watts of Power.
- UV Lighting Strip ensures VHP & Formalin Sterilisation.
- Electric Sash Adjustment & Closure.
- H14 HEPA Filters that help create an ISO Class 4 Clean Environment.
- Air Diffuser and Sloped Intake for efficient and smooth airflow.
- Mirror Finish 316 Grade Stainless Steel Work Surface.
- Infra-Red Intelligent Sash Safety (IRISS).

Visionaire[®] Touchscreen Control

Our revolutionary Visionaire[®] Touchscreen Control is designed to control & monitor all aspects of product operation including; a Visual Display of Inlet and Downward Air Velocity (m/sec or ft/min), Sash Adjustment, Control for: UV Light, Electrical Socket, Cabinet Power and Light Variable intensity.

Audible and Visual alarms alert the user to maintenance requirements including: Service, UV Light replacement and Cabinet Hours Run. Formalin and VHP Cycles are also controllable.

Size Specification (W x D x H)

MSC 1200 C1

1200mm x 700mm x 1320mm