

ELIBATH



TISSUE WATERBATH



Thanks to its elongated shape, EliBath can be oriented in any direction. This helps the technician to move freely and naturally, and improves the work quality and efficiency.

Why choose EliBath?

***Ergonomic shape. Small footprint.**

Minimizes the space occupied on the workbench, while maintaining a large working area

***Digital control.**

Accurate operation and automatic water temperature control

***Easy water replacement.**

Lightweight, easily transported and handled

***Adjustable light setting.**

9 different LED colors, adjustable brightness





Optimize your working station!



EliBath is the new concept of the tissue waterbath. The elongated shape allows to work across a preferred orientation for picking up sections. A full usable depth of 65 mm makes it possible to insert slides vertically. Thanks to its innovative design, EliBath maximizes the working volume, while being extremely compact: 1,8 liter of water and a useful depth of 65 mm for just 1.4 Kg of weight and a footprint of 29 x 20 cm. The magnetic power port allows fast connection and disconnection of the unit, which, together with its lightweight, allows water to be replaced often, avoiding contaminations risks. The waterbath can be controlled by a large external display, which also allows to adjust the light settings and to activate advanced functions.

Technical specifications

Temperature control

Temperature range	RT to 70 C
Voltage	24 V DC
Power	Max 280 W
Capacity	1.8 L

Weight

Waterbath	1.4 Kg
-----------	--------

Dimensions

Total:	290 x 200 x 87 mm
Tank:	243 x 142 x 67 mm
Max usable depth:	65 mm
Filling level:	60 mm

Connectivity

LabIOT system compatible

